Tailored Implants made of PVDF

DynaMesh[®]

- (de) Wiederverwendbare Mehrweginstrumente aus chirurgischem Stahl
- (en) Reusable instruments made of surgical steel
- (cs) Opakovaně použitelné nástroje z chirurgické oceli
- (el) Επαναχρησιμοποιήσιμα εργαλεία πολλαπλών χρήσεων από χειρουργικό χάλυβα
- (es) Instrumentos reutilizables de acero quirúrgico
- (fr) Instruments réutilisables en acier chirurgical
- (it) Strumenti riutilizzabili in acciaio chirurgico
- (nl) Herbruikbare instrumenten van chirurgisch staal
- (no) Gjenbrukbare instrumenter av kirurgisk stål
- (pt) Instrumentos reutilizáveis de aço cirúrgico
- (ro) Instrumente reutilizabile din oțel chirurgical
- (zh-cn) 手术钢制成的可重复使用器械
- (zh-tw) 手術鋼製成的可重複使用的器械

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DynaMesh®-ISR01
DynaMesh®-IVT01 / -IVT02
DynaMesh®-IST01 / -IST02 / -IST03

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DynaMesh®-ISR01
DynaMesh®-IVT01 / -IVT02
DynaMesh®-IST01 / -IST02 / -IST03

1. DESCRIPTION

DynaMesh® instruments are reusable instruments made of surgical steel.

2. FIELDS OF APPLICATION

The **DynaMesh®** instruments are used together with the urological/gynaecological products from the **DynaMesh®** family for the surgical treatment of urinary incontinence and in prolapse surgery.

3. WARNING NOTICES

When using the instruments, punctures or tears to the organs, blood vessels or nerves may occur, which then require surgical treatment.

4. CLEANING AND STERILISATION

DynaMesh® instruments are supplied unsterile. The instruments must be instrumentally and/or manually cleaned, checked for operability and sterilised prior to their first use and further applications.

RECOMMENDATIONS FOR STERILISATION:

Vapour sterilisation: Any combination of time and temperature providing a Sterility Assurance Level (SAL) of at least 10^{-6} , e.g. the recommendation according to the European Pharmacopoeia (121°C, 15 min)

EO sterilisation: Any method validated to a SAL of at least 10^{-6} .

Pursuant to the recommendations of the Robert Koch Institute, the instruments are to be discarded in case of a suspected contamination with prions.

5. STORAGE CONDITIONS

DynaMesh® must be stored in a dry and clean state in appropriate packages.

6. EMBODIMENTS

DynaMesh® instruments are available in the following versions:

DynaMesh® -ISR01: Instrument for access for retrosymphyseal tape insertion and for transvaginal access for transgluteal mesh insertion.

DynaMesh® -IVT01 / 02: Instrument for transvaginal access for transobturatoric mesh insertion, for access for transobturatoric tape insertion and for laparotomic access for retroperitoneal tape insertion.

DynaMesh® -IST01 / 02 / 03: Instrument set comprising 2 instruments (right- and left-hand side) for transvaginal access for transobturatoric mesh insertion and for access for transobturatoric tape insertion.

DynaMesh® instruments comply with the basic requirements of the directive 93/42/EWG concerning medical devices.

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Products:

All reusable surgical products manufactured
by FEG

Warning

Remove and dispose of tip guard before first use!
Clean instruments paying particular attention to
the lumens. Use appropriate brushes.

Limitation of processing

Frequent reprocessing has slight effects on the
instrument. The end of the product life usually
depends on the wear and damage due to use.

User/point of use

Remove surface soiling with a disposable cloth.
Store instruments properly for disposal. Do not
overload the disposal tray.

Containment and transport

No special requirements.
It is recommended that instruments are
reprocessed immediately after use. Dry
disposal in a closed system is preferable.

Preparation for cleaning

No dismantling involved.
In case of heavy soiling, the instruments can
safely be pre-cleaned in an ultrasonic bath
using a VAH-listed detergent/disinfectant.

Automated cleaning

Equipment: Miele washer and disinfector, type
7828

Detergent: Sekumatik MultiClean 0.6%,
Ecolab, mildly alkaline

Neutraliser: Sekumatik FNP 0.08%, Ecolab

Loading:

Place instrument in a tray and position side by
side in the washer basket
Take care not to overload the trays or block the
spray arms.

Remove the instruments immediately after the
end of program to avoid corrosion caused by
residual moisture in the closed washer.

Vario program cycle:

Pre-wash (3 min, cold), cleaning (10 min, warm

55°C), neutralising (2 min, warm), rinsing (1 min,
cold, demineralised water), thermal disinfection
(5 min, 93°C), drying (30 min, 115°C)

Manual cleaning

Equipment: Ultrasonic bath, cleaning and
disinfecting chemical: Sekusept aktiv 2%,
Ecolab, brush, demineralised water

- Place instruments in the cleaning and
disinfecting solution.
- Treat with ultrasound for 5 minutes.
- Use a soft plastic brush to clean the
instrument (lumen).
- Leave instrument to soak in the solution for
15 min as per instructions provided by the
process chemical manufacturer.
- Rinse thoroughly with plenty of
demineralised water to remove any
remaining dirt particles.
- Dry the instruments immediately after
cleaning.

Care

No care products required.

Cleanliness check and functional test

The instruments must be macroscopically
clean, i.e. free from visible residue. Inspect
the instruments carefully, paying particular
attention to critical areas such as the lumens.
Damaged and worn instruments must be
replaced.

Packaging

Individually or in the tray. Packaging material:
set/container, fleece, laminate foil pack

Sterilisation

Steam sterilisation at 134°C, 18 min

Storage

Dry and clean in suitable packaging

The instructions provided above have been validated
by the manufacturer as being CAPABLE of preparing
the product for reuse. It remains the responsibility
of the processor to ensure that the reprocessing
as actually performed using equipment, materials
and personnel in the reprocessing facility achieve
the desired result. This normally requires validation
and routine monitoring of the process. Furthermore,
the product should be carefully inspected by the
processor for deviations from the instructions
provided to evaluate their effectiveness and potential
detrimental effects.

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