
Expert Technologies in PVDF

# DynaMesh<sup>®</sup> - CESA/-CERESA

# DynaMesh<sup>®</sup> - VASA/-VARESA

de	Implantate für die Beckenbodenchirurgie
en	Implants for pelvic floor surgery
ar	غرسات من أجل جراحة قاعدة الحوض
cs	Implantáty pro chirurgii pánevního dna
da	Implantater til bækkenbundskirurgien
es	Implantes para cirugía de suelo pélvico
fa	ایمپلنت برای جراحی کف لگن
fi	Istukkeet lantionpohjan kirurgiaan
fr	Implants pour la chirurgie du plancher pelvien
it	Impianti per la chirurgia del pavimento pelvico
he	שתל לכירורגיה של רצפת האגן
ko	골반저 수술용 임플란트
nl	Implantaten voor de bekkenbodemchirurgie
no	Implantater for bekkenbunnskirurgi
pl	Wszczepy do chirurgii dna miednicy
pt	Implantes para cirurgia do solo pélvico
ro	Implanturi pentru planșeul pelvin
sk	Implantáty pre chirurgiu panvového dna
sv	Implantat för bäckenbottenkirurgi
tr	Pelvik taban cerrahisi için implantlar
zh	盆底手术植入补片

	①	②	③	④
de	Einführungshilfen	Stichmarkierungen	Fixationsbereiche Rektopexie	Fixationsbereich Vaginal-/Cervikal-Stumpf
en	Insertion aids	Stitch marks	Fixation areas for rectopexy	Fixation area for vaginal/cervical stump
ar	وسائل الإدخال المساعدة	علامات الغرزات	مناطق تثبيت المستقيم	منطقة تثبيت رباط الرحم / رباط عنق الرحم
cs	Pomocné zaváděcí části	Stehové značky	Fixační oblasti rektopexie	Fixační oblast poševního/cervikálního pahýlu
da	Indføringshjælp	Stikmarkeringer	Fikseringsområder rektopexi	Fikseringsområde vaginal/cervikal-stump
es	Ayudas de inserción	Marcas de punción	Zonas de fijación para retopexia	Zona de fijación para muñón vaginal / cervical
fa	نوار ورودی	محل بخیه	محل تثبیت رکتوپکسی	محل تثبیت استامپ واژن و یا سرویکال
fi	Sisäänvientiöjaimet	Pistomerkinnät	Kiinnitysalueet rektopeksiassa	Kiinnitysalue emätin/kohdunkaulatyngässä
fr	Dispositifs d'insertion	Repères de suture	Zones de fixation retopexie	Zone de fixation Dôme vaginal / Moignon cervical
it	Ausili di introduzione	Contrassegni	Area di fissaggio per retopessi	Area di fissaggio per calotta vaginale/moncone cervicale
he	סיוע בהשתלה	סימוני תפירה	אזור הקיבוע ברקטופוקסיה	אזור הקיבוע לציחת נרתיק/ צוואר הרחם
ko	삽입 보조 도구	스티치 표시	직장 고정술 고정 영역	자궁 경부/질 처단부 고정 영역
nl	Inbreng hulpmiddelen	Steekmarkeringen	Fixatiebereiken rektopexie	Fixatiebereik vaginale apex/cervixstomp
no	Innføringshjelpere	Stingmarkeringer	Fikseringsområder rektopeksi	Fikseringsområde vaginal/cervix-stump
pl	Taśmy wprowadzające	Punkty założenia szwów	Obszary mocowania rektopeksji	Obszar mocowania kikuta pochwy/szyjki macicy
pt	Auxiliares de inserção	Marcações de pontos	Áreas de fixação para retopexia	Área de fixação para colo do útero/cúpula vaginal
ro	Dispozitive ajutoare de introducere	Marcaje de sutură	Zone de fixare retopexie	Zone de fixare apex vaginal/de col uterin
sk	Zavádzacie pomôcky	Značky stehov	Oblasti fixácie pri rektopexii	Oblasť fixácie na vaginálny pahýl / pahýl krčku maternice
sv	Införingshjälpmedel	Stickmarkeringar	Fixeringsområde rektopexi	Fixeringsområde slidbotten/livmodertapp
tr	Kılavuzlar	Dikiş işaretleri	Rektopeksi için sabitleme bölgeleri	Vajinal/servikal güdük sabitleme bölgesi
zh	导引辅助端	缝线标记	直肠固定术固定区域	阴道/宫颈残端固定区域

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## Instructions for use

### 1. Description

**DynaMesh® CESA/CERESA** and **VASA/VARESA** are used to reinforce connective tissue structures and ligaments. They do not have any pharmacological effect of their own and are non-toxic.

The implants are knitted from non-absorbable, biostable polyvinylidene fluoride monofilament. For handling and recognition purposes, some of the threads are coloured black. These threads also make it possible to visualise the mesh postoperatively by MRI.

The specially developed knitted structure and the optimal pore size of the implants are perfectly matched to the fields of application.

### 2. Fields of application

The **DynaMesh® CESA/CERESA** and **VASA/VARESA** implants have been specially developed for pelvic floor reconstruction, and particularly for reinforcing or replacing the uterosacral ligaments. The implants are used in the treatment of a prolapse of the internal genitalia, such as prolapse of the vaginal stump or uterus, and the resulting urge incontinence.

The surgical technique represents a modified version of the abdominal sacropexy, in which the uterosacral ligaments are reinforced or replaced by the implant on both sides. Different versions of **DynaMesh®** implants are available according to the type of operation being performed and the anatomical defect in the patient:

**DynaMesh® CESA** is used to affix the cervical stump to the sacrum  
**DynaMesh® CERESA** is used to affix the cervical stump to the sacrum and for additional rectopexy in cases of simultaneous faecal incontinence

**DynaMesh® VASA** is used to affix the vaginal stump to the sacrum  
**DynaMesh® VARESA** is used to affix the vaginal stump to the sacrum and for additional rectopexy in cases of simultaneous faecal incontinence

For the exact dimensions and structures of **DynaMesh®** implants, please refer to the current product range.

### Fixation

The implants are affixed to the sacrum on both sides at the stitch marks. On the vaginal or cervical stump, the mesh should be fixed at four points, as smoothly as possible. In both cases we recommend using non-absorbable suture material. In the case of rectopexy, the mesorectum should be affixed to the fixation part of the implant with absorbable suture material. We have no experience of using other fixation materials such as staplers or glues.

### 3. Contraindications

All the usual contraindications relating to the use of implants for tissue support also apply to **DynaMesh® CESA/CERESA** and **VASA/VARESA** implants.

**DynaMesh®** implants should not be used in patients with obvious infections.

### 4. Warnings

**DynaMesh® CESA/CERESA** and **VASA/VARESA** implants are intended for single use only. Re-use can lead to functional, immunological, toxicological or hygienic impairments. Excess mesh on the vaginal/cervical stump can be cut off after suturing. When doing so, please take care not to damage the three central load-bearing structures, which are coloured black. The insertion aids must be cut off approximately 0.5 cm behind the stitch marks towards the end of the band. They must not remain in the body.

The user must be familiar with the surgical technique involved in the placement of mesh before using **DynaMesh® CESA/CERESA** and **VASA/VARESA**.

If **DynaMesh®** implants are used in a contaminated wound, a subsequent infection is possible which may necessitate removal of the material.

Treatment with the **DynaMesh® CESA/CERESA** and **VASA/VARESA** may bring to light a masked stress incontinence, which might require additional treatment.

The known risks in prolapse surgery with the aid of alloplastic materials include erosions, infections, pain, bleeding, dyspareunia, organ perforation, voiding dysfunction, de novo stress incontinence and shrinkage of the mesh.

The following complications may arise with alloplastic materials in general: seromas, haematomas, fistula formation, inflammations, recurrences or arrosions.

Alloplastic material should only be used after carefully considering conventional surgical procedures and fully informing the patient. The patient should be advised that physical activities such as heavy lifting, jogging, riding or sexual intercourse must only be resumed after consultation with the doctor in attendance.

Strictly sterile handling must be ensured up to implantation. The procedure must only be carried out in aseptic rooms.

Damaged implants must no longer be used. Contact with sharp objects in particular must be avoided.

**DynaMesh® CESA/CERESA** and **VASA/VARESA** must not be implanted if:

- the expiry date has passed or
- the packaging is damaged or
- there is reason to suspect the packaging is not intact.

The implant is supplied in a sterile double packaging. The inner packaging should not be opened until the mesh is about to be used. It is recommended that the mesh be handled only with sterile, unpowdered gloves and clean instruments.

The surgeon is responsible for complications arising from an incorrect indication, an inadequate surgical technique or failure to ensure sterile conditions. In all cases, it is recommended that the operation be followed up with appropriate post-operative treatment.

### 5. Side effects

There are no indications in the literature that the polyvinylidene fluoride used has any adverse side effects. As with all foreign bodies, **DynaMesh®** implants may have a negative effect on a pre-existing infection and may lead to temporary local irritations.

### 6. Sterility

**DynaMesh®** implants are sterilised with ethylene oxide.

**DynaMesh®** must never be resterilised.

The products must not be used if the packaging has been damaged or opened!

Opened, unused products must be discarded!

### 7. Storage conditions

The recommended storage conditions are:

- Store in a dry place at room temperature.
- Protect from moisture and the direct effect of heat.
- Store only in the original packaging.

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