Tailored Implants made of PVDF

DynaMesh® - LICHTENSTEIN

- (de) Netzimplantat für die Hernienchirurgie aus PVDF-Monofilament
- (en) Mesh implant for hernia surgery, made of PVDF monofilament
- (ar) زرع شبكي خاص بجراحة الفتق من خيط وحيد مكون من ثنائي الفلوريد المتعدد الفينيليدي
- (cs) Síťkový implantát pro herniální chirurgii z PVDF-monofilamentu
- (da) Netimplantat af PVDF-monofilament til herniekirurgi
- (el) Εμφύτευμα πλέγματος για χειρουργική επέμβαση κήλης, κατασκευασμένο από μονονημάτιο PVDF
- (es) Implante de malla para cirugía de hernias fabricado en monofilamento de PVDF
- (fi) PVDF-monofilamentista valmistettu verkkoimplantti tyräkirurgiaan
- (fr) Filet prothétique pour la chirurgie des hernies, en monofil de PVDF
- (it) Reti per la chirurgia dell'ernia in monofilamento di PVDF
- (nl) Netimplantaat voor de herniachirurgie uit PVDF-monofilament
- (no) Nettingimplantat i PVDF-monofilament for brokkkirugi
- (pl) Siatka chirurgiczna z włókien monofilamentowych z polifluorku winylidenu (PVDF) do chirurgicznego leczenia przepuklin
- (pt) Implante de rede de monofilamentos de PVDF para a cirurgia de hérnias
- (ro) Implant de plasă din monofilament PVDF ptr.chirurgia herniei
- (ru) Сетчатый имплантат для герниопластики из ПВДФ монопнитей
- (sk) Nettingimplantat i PVDF-monofilament for brokkkirugi
- (sv) Nätimplantat för bräckkirurgi av PVDF-monofilament
- (tr) Herni cerrahisinde kullanılan PVDF Monofilament ağ implant
- (zh) 疝氣手術植入網片 (由 PVDF 單纖維製成)

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Instructions for Use

1. Description

DynaMesh®-LICHTENSTEIN is a mesh implant for reinforcing connective tissue structures. **DynaMesh®-LICHTENSTEIN** does not have any pharmacological effect of its own and is non-toxic. **DynaMesh®-LICHTENSTEIN** is knitted using non-absorbable, biostable polyvinylidene fluoride monofilament. For handling and recognition purposes, the threads are partially coloured green or black. The coloured marking threads are used to improve placement of the mesh during the operation, in addition, the black threads enable postoperative recognition of the mesh in an MRI (Magnetic Resonance Imaging) scan. The specially developed knitted structures and optimal mesh size of **DynaMesh®-LICHTENSTEIN** are perfectly matched to the dynamometry of the abdominal wall, in particular the inguinal region.

2. Fields of application

DynaMesh®-LICHTENSTEIN serves to support the tissue and stabilise the fascial structures of the inguinal region. **DynaMesh®-LICHTENSTEIN** was specially developed for the conventional repair of inguinal hernias according to Lichtenstein. **DynaMesh®-LICHTENSTEIN** must not be used in intraperitoneal onlay mesh technique (IPOM technique).

DynaMesh®-LICHTENSTEIN comes in a size which fits the majority of patients and can therefore be used immediately without needing to be cut to size. The mesh can be cut to size to cater for anatomical peculiarities in individual cases of course, but this is not generally recommended. The dimensions correspond to those in the original description of the procedure according to Lichtenstein. Should it be necessary to adjust the size in rare cases, the mesh can be cut with sharp scissors or a scalpel outside the operating field.

Fixation

It is recommended that the mesh is fixed according to Lichtenstein's original method, using a continuous suture on the inguinal ligament and interrupted sutures in a cranial and lateral direction. The suture thickness is at the surgeon's discretion. Alternatively, clips, staples or biological or synthetic adhesives can be used, provided appropriate precautions are taken. This differs from Lichtenstein's original method however.

3. Contraindications

All the usual contraindications relating to the use of implants for tissue support also apply to the application of **DynaMesh®-LICHTENSTEIN**. Among others, these include:

- The mesh should not be used where there has been a previous enterotomy, due to the increased risk of infection.
- The mesh should be separated from the abdominal cavity by the peritoneum.
- The mesh is not suitable for use as a plug.
- An inguinal hernia repair using a mesh is totally inappropriate where children are concerned and caution at least must be exercised when using a mesh in patients who have obvious infections or who are pregnant.

4. Warnings

DynaMesh®-LICHTENSTEIN is intended for single use only. Re-use can lead to functional, immunological, toxicological or hygienic impairments. Users must be familiar with the technique

used to repair inguinal hernias before applying **DynaMesh®-LICHTENSTEIN**. If **DynaMesh®-LICHTENSTEIN** is used in a contaminated wound, a subsequent infection is possible which may make it necessary to remove the material. The implantation must be carried out under strictly sterile conditions. The application must be carried out only in aseptic rooms. Damaged implants must no longer be used. Contact with sharp objects in particular must be avoided.

DynaMesh®-LICHTENSTEIN must not be implanted if:

- the sell-by date has expired,
- the packaging is damaged or
- there is reason to suspect the packaging is not intact.

The implant is supplied in sterile double packaging. The inner packaging should not be opened until the mesh is about to be used. It is recommended that the mesh be handled only with sterile, unpowdered gloves and clean instruments.

The following complications may arise with meshes in general: seromas, haematomas, relapses and infections, as well as visceral adhesions, should the implant come into direct contact with the intestines. The surgeon is responsible for complications arising from an incorrect indication, an inadequate operating technique or failure to ensure sterile conditions. It is recommended in all cases that the operation be followed up with appropriate post-operative treatment.

5. Side effects

There are no indications in the literature that the polyvinylidene fluoride used has adverse side effects. Like all foreign bodies, **DynaMesh®-LICHTENSTEIN** may have a negative effect on a pre-existing infection or lead to temporary local irritations.

6. Sterility

The product is sterilised with ethylene oxide. **DynaMesh®-LICHTENSTEIN** must on no account be resterilised. The product must not be used if the packaging has been damaged or opened. Opened, unused products must be discarded.

7. Storage conditions

The recommended storage conditions are:

- Store in a dry place at room temperature.
- Protect from moisture and the direct effect of heat.
- Store only in the original packaging.
- Do not use after the expiry date.

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