Tailored Implants made of PVDF

DynaMesh[®] -PR / soft / visible
DynaMesh[®] -PRR soft / visible
DynaMesh[®] -PRS soft / visible
DynaMesh[®] -PR2 soft / visible
DynaMesh[®] -PR4 soft / visible

de	Netzimplantate für die Beckenbodenchirurgie
en	Mesh implants for pelvic floor surgery
cs	Sít'kové implantáty pro chirurgii pánevního dna
el	Εμφυτεύματα πλέγματος για χειρουργικές επεμβάσεις του πυελικού εδάφους
es	Implantes de malla para la cirugía del suelo pélvico
it	Impianti a rete per la chirurgia del pavimento pelvico
nl	Meshimplantaten voor de bekkenbodenchirurgie
pt	Implantes têxteis para cirurgia do solo pélvico
ro	Implant de plasă din monofilament PVDF ptr.chirurgia herniei
tr	Pelvik taban cerrahisi için ağ implantlar
zh	盆底手術網狀植入補片

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DynaMesh®-PR / soft / visible
DynaMesh®-PRR soft / visible
DynaMesh®-PRS soft / visible
DynaMesh®-PR2 soft / visible
DynaMesh®-PR4 soft / visible

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1. DESCRIPTION

DynaMesh®-PR are mesh implants for reinforcing connective tissue structures and ligaments. **DynaMesh®-PR** mesh implants do not have any pharmacological effect of their own and are non-toxic. The implants are knitted using non-absorbable, biostable polyvinylidene fluoride monofilaments. For handling and recognition purposes, the threads are partially coloured green or black. The coloured marking threads are used to improve placement of the meshes during the operation. In addition, the black threads enable postoperative recognition of the mesh in an MRI scan. The specially developed knitted structure and the optimal mesh size of **DynaMesh®-PR** mesh implants are perfectly matched to the fields of application.

2. FIELDS OF APPLICATION

DynaMesh®-PR mesh implants have been specially developed for pelvic floor reconstruction with mesh implants. The implants serve to support and stabilise fascial structures, connective tissue and ligaments. **DynaMesh®-PR** mesh implants are used, for example, in the surgical treatment of prolapse of the vaginal stump or uterine prolapse, as well as in the treatment of a rectocele, enterocele or cystocele.

The choice of the operation technique and the final position of the mesh is left to the surgeon's discretion. Depending on the type of operation being carried out and the anatomical defect in the patient, it is possible to choose between various versions of **DynaMesh®-PR** mesh implants:

DynaMesh®-PR / soft / visible, DynaMesh®-PRR soft / visible

Prolapse of the vaginal stump or uterine prolapse:

Sacrocolpopexy/open or laparoscopic technique

DynaMesh®-PRS soft / visible

Prolapse of the vaginal stump: Sacrocolpopexy/open or laparoscopic technique

DynaMesh®-PR2 soft / visible

Rectocele, enterocele: Transvaginal mesh plasty (posterior)

DynaMesh®-PR4 soft / visible

Cystocele: Transvaginal mesh plasty (anterior)

For the exact dimensions and structures of **DynaMesh®-PR** mesh implants, please refer to the current product range. When using **DynaMesh®-PR** mesh implants, it is possible to choose between various versions of insertion instruments depending on the approach and the position of the mesh:

DynaMesh®-ISR01: Instrument for transvaginal access for retrosymphyseal tape insertion and for transvaginal access for transgluteal mesh insertion.

DynaMesh®-IVT01: Instrument for transvaginal access for transobturatoric mesh insertion, for transvaginal access for transobturatoric tape insertion and for laparotomic access for retroperitoneal tape insertion.

DynaMesh®-IST01 / 02 / 03: Instrument set comprising 2 instruments (right- and left-hand side) for transvaginal access for transobturatoric mesh insertion and for transvaginal access for transobturatoric tape insertion.

3. CONTRAINDICATIONS

All of the usual contraindications relating to use as an implant for tissue support also apply to the use of **DynaMesh®-PR** mesh implants.

DynaMesh®-PR mesh implants should not be used in patients with obvious infections, if the patient is currently or intends to become pregnant, or if the patient is still growing.

4. WARNINGS

DynaMesh®-PR mesh implants are intended for single use only. Re-use can lead to functional, immunological, toxicological or hygienic impairments. The user must be familiar with the surgical technique involved in the placement of mesh before using **DynaMesh®-PR**. The use of alloplastic mesh should only be undertaken after a careful consideration of classic surgical procedures has been made and the patient has been offered a comprehensive explanation. If **DynaMesh®-PR** mesh implants are used in a contaminated wound, a subsequent infection is possible which may make it necessary to remove the material. The known risks in descensus and prolapse surgery with the aid of alloplastic materials include erosions, infections, pain, bleeding, dyspareunia, organ perforation, difficulty urinating and shrinkage of the mesh. The following complications may arise with meshes in general: Seromas, haematomas, fistula formations, inflammations, recurrences or arrosions. The patient should be advised that physical activities such as heavy lifting, jogging, riding or sexual intercourse should be resumed only after consultation with the doctor in attendance. The procedure must be carried out under strictly sterile conditions up to implantation. The application must be carried out only in aseptic rooms. Damaged implants must no longer be used. Contact with sharp objects in particular must be avoided.

DynaMesh®-PR mesh implants must not be implanted if,

- the expiry date has passed or
- the packaging is damaged or
- there is reason to suspect the packaging is not intact.

The implant is packaged in a sterile double packaging. The inner packaging should not be opened until the mesh is about to be used. It is recommended that the mesh be handled only with sterile, unpowdered gloves and clean instruments. The surgeon is responsible for complications arising from an incorrect indication, an inadequate operating technique or failure to ensure sterile conditions. In all cases, it is recommended that the operation be followed up with appropriate post-operative treatment.

5. SIDE EFFECTS

There are no indications in literature that the polyvinylidene fluoride used has adverse side effects. As with all foreign bodies, **DynaMesh®-PR** mesh implants may have a negative effect on a pre-existing infection and may lead to temporary local irritations.

6. STERILITY

DynaMesh®-PR mesh implants are sterilised with ethylene oxide. **DynaMesh®-PR** must on no account be resterilised! The products must not be used if the packaging has been damaged or opened!

Opened, unused products must be discarded!

7. STORAGE CONDITIONS

The recommended storage conditions are:

- Store in a dry place at room temperature.
- Protect from moisture and the direct effect of heat.
- Store only in the original packaging.

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