Tailored Implants made of PVDF

# DynaMesh®-PRM

de	Netzimplantate für die Beckenbodenchirurgie
en	Mesh implants for pelvic floor surgery
cs	Sít'kové implantáty pro chirurgii pánevního dna
el	Εμφυτεύματα πλέγματος για χειρουργικές επεμβάσεις του πυελικού εδάφους
es	Implante de malla para la cirugía del suelo pélvico
fr	Implants textiles en filet pour la chirurgie du plancher pelvien
it	Impianti di rete per la chirurgia del pavimento pelvico
nl	Meshimplantaten voor de bekkenbodenchirurgie
no	Nettingimplantater for bekkenbunnskirurgi
pt	Implantes têxteis para cirurgia do solo pélvico
ro	Implant de plasă ptr. chirurgia peretelui pelvic
zh-tw	用於盆腔底壁手術的植入網片

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- DynaMesh®-PRM
- DynaMesh®-PRM visible

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- DynaMesh®-PRM
- DynaMesh®-PRM visible

## 1. DESCRIPTION

**DynaMesh®-PRM** is a mesh implant for reinforcing connective-tissue structures and ligaments, in particular for the treatment of male urinary stress incontinence.

**DynaMesh®-PRM** does not have any pharmacological effect of its own and is non-toxic.

**DynaMesh®-PRM** is knitted from non-absorbable, bio-stable polyvinylidene fluoride monofilament fibre. For handling and recognition purposes, the threads are partially coloured green or black. The black coloured thread material in the mesh structure enables post-operative recognition of the mesh in an MRI (Magnetic Resonance Imaging) scan. The specially developed knitted structures and the ideal mesh size of **DynaMesh®-PRM** are optimally matched to the fields of application for which it is intended.

## 2. FIELDS OF APPLICATION

**DynaMesh®-PRM** is used to support and stabilise connective-tissue structures and ligaments. Common applications are suburethral slings for treating male urinary stress incontinence.

The surgeon selects the method of surgery and position of the sling, both of which will depend on his/her preference, on the one hand, and on the instruments used, on the other. These instruments may be procured separately. Normally, the mesh implant is implanted through a perineal access. The bulbospongiosus muscle is divided in the midline and the urethral bulb is detached from the fibrous nucleus until the urethral bulb is mobile. The tape is inserted transobturarily with a tunneler. The mesh part is fixed to the mobile part of the urethral bulb using sutures. Traction on the tape draws the bulb in a cranial direction in the retrourethral space, which should correct suburethral relaxation to strengthen the sphincter.

When the tape is finally in position, a tunnel is formed with an Overholt from the perineal wound to the incision site. The ends of the tape are pulled back and cut off at the level of the wound edges. To prevent the tape being displaced in the first weeks of wound healing, fixation of the ends of the tape is necessary. For this purpose, the ends of the tape are fixed with a synthetic adhesive near the exit point of the ends of the tape. Alternatively, the ends of the tape can be loosely connected with a long-term, absorbable suture material with a thickness of 2-0 at the incision site level.

The user must be familiar with the technique involved in implanting a suburethral sling before fitting a **DynaMesh®-PRM**.

## 3. CONTRAINDICATIONS

All the usual contraindications for use as a tissue support implant also apply to the use of **DynaMesh®-PRM**. In particular the implant should not be used where urinary tract infections are present or in patients with blood coagulation disorders.

## 4. WARNINGS

**DynaMesh®-PRM** is intended for single use only. Re-use can lead to functional, immunological, toxicological or hygienic impairments.

If **DynaMesh®-PRM** is used in a contaminated wound, there is a possibility of subsequent infection, which may necessitate removal of the material.

Insertion of the needle may cause damage to blood vessels, nerves, the bladder or intestine, which will require immediate surgical attention.

Meshes in general may lead to the following complications: seromas, haematomas, infections, fistulisation, inflammations or arrosions.

Up to implantation strict care must be taken to ensure sterile handling. Application may only take place in aseptic rooms.

Damaged implants must not be used; in particular, contact with sharp objects must be avoided.

**DynaMesh®-PRM** must not be implanted if

- the minimum shelf life has been exceeded or
- the packaging is damaged or
- there is any doubt about its integrity.

The implant is packed in sterile double-packaging. The inner packaging should not be opened until the mesh is about to be used. It is recommended that the mesh only be handled with sterile, unpowdered gloves and clean instruments.

The reusable instruments are to be checked before use to ensure that they are in perfect condition.

The surgeon is responsible for any complications which may be due to an incorrect indication, an inadequate operating technique or an error in maintaining aseptic conditions. In all cases it is recommended that the intervention is accompanied by appropriate follow-up.

The patient's attention should be drawn to the fact that such physical activities as heavy lifting, jogging, riding or sexual intercourse should not be resumed until the attending physician has been consulted.

## 5. SIDE-EFFECTS

There are no references in the literature to any undesirable side-effects caused by the polyvinylidene fluoride used. As with all foreign bodies, **DynaMesh®-PRM** can have a detrimental effect on a pre-existing infection and lead to temporary local irritation.

## 6. STERILITY

The product is sterilised with ethylene oxide.

**DynaMesh®-PRM** must not be re-sterilised!

The product must not be used if the packaging is damaged or already opened!

Any products which have been opened but not used must be discarded!

## 7. STORAGE CONDITIONS

The recommended storage conditions are as follows:

- Store in a dry place at room temperature
- Protect against moisture and direct heat
- Only store in the original packaging

Do not use if the expiry date has elapsed!

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