
DynaMesh[®] - PRP soft / visible

de	Implantate für die Beckenbodenchirurgie	3
en	Implants for pelvic floor surgery	4
br	Импланти за хирургията на тазовото дъно	5
cs	Implantáty pro chirurgii pánevního dna	6
da	Implantater til bækkenbundskirurgien	7
el	Εμφυτεύματα για χειρουργικές επεμβάσεις του πυελικού εδάφους	8
es	Implantes para cirugía de suelo pélvico	9
et	Vaagnarõhjakirurgia implantaadid	10
fi	Istukkeet lantionpohjan kirurgiaan	11
fr	Implants pour la chirurgie du plancher pelvien	12
hr	Implantati za kirurgiju dna zdjelice	13
hu	Implantátumok medencealapi sebészeti beavatkozásokhoz	14
it	Impianti per la chirurgia del pavimento pelvico	15
lt	Implantai dubens dugno operacijoms	16
lv	Implanti iegurņa pamatnes ķirurģijai	17
nl	Implantaten voor de bekkenbodemchirurgie	18
no	Implantater for bekkenbunnskirurgi	19
pl	Implanty do chirurgii dna miednicy	20
pt	Implantes para cirurgia do solo pélvico	21
ro	Implanturi pentru planșeul pelvin	22
ru	Имплантаты для хирургии тазового dna	23
sl	Vsadki za medenično dno	24
sr	Implantati za hirurgiju karličnog dna	25
sv	Implantat för bäckenbottenkirurgi	26
tr	Pelvik taban cerrahisi için implantlar	27

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- **DynaMesh®-PRP soft / visible** 03 cm x 15 cm
- **DynaMesh®-PRP visible** 03 cm x 18 cm
- **DynaMesh®-PRP visible** 17 cm x 15 cm

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1. DESCRIPTION

DynaMesh®-PRP soft / visible are mesh implants for reinforcing connective tissue structures and ligaments. **DynaMesh®-PRP soft / visible** mesh implants do not have any pharmacological effect of their own and are non-toxic. The implants are knitted from non-absorbable, biostable polyvinylidene fluoride monofilament. For handling and recognition purposes, the threads are partially coloured green or black. The coloured marking threads are used to improve placement of the meshes during the operation. In addition, the black threads enable postoperative recognition of the mesh in an MRI scan. The specially developed knitted structure and the optimal mesh size of **DynaMesh®-PRP soft / visible** mesh implants are perfectly matched to the fields of application.

2. FIELDS OF APPLICATION

DynaMesh®-PRP soft / visible mesh implants have been specially developed for pelvic floor reconstruction with mesh implants. The implants serve to support and stabilise fascial structures, connective tissue and ligaments.

DynaMesh®-PRP soft / visible mesh implants are used in the surgical treatment of a prolapse of the vaginal stump or uterine prolapse in the pectopexy technique. With an accompanying prolapse of the anterior and / or posterior vaginal wall (pulsion cystocele / rectocele), additional stabilisation of the affected vaginal wall can be achieved with the **DynaMesh®-PRP visible** implants measuring 17 cm x 15 cm.

Following hysterectomy:

To this end, the iliopectineal ligament is dissected on both sides in the area between the vena iliaca externa and the medial umbilical ligament, and the peritoneum is opened up to the vaginal or cervical stump. The middle **DynaMesh®-PRP soft / visible** mesh segment is fixed to the vaginal or cervical stump, the implant ends are fixed to the lateral part of the iliopectineal ligament (where the psoas muscle crosses / level with first sacral vertebra). The peritoneum is then completely closed.

Two implant sizes are available in the dimensions **DynaMesh®-PRP soft / visible** 03 cm x 15 cm and **DynaMesh®-PRP visible** 03 cm x 18 cm. In the majority of cases, **DynaMesh®-PRP soft / visible** 03 cm x 15 cm is used. With greatly shortened vaginas, e.g., following a radical hysterectomy, **DynaMesh®-PRP visible** 03 cm x

18 cm can be optionally used.

In the case of uterine preservation: Proceed in the same manner, with fixation to the cervix. With a normal sized uterus, the **DynaMesh®-PRP visible** implants measuring 03 cm x 18 cm should be used and fixed in place on the posterior cervix. With smaller uteri (below 100 g), anterior fixation of the **DynaMesh®-PRP soft / visible** implants measuring 03 cm x 15 cm can be selected as an alternative.

With accompanying cystocele or rectocele: With an accompanying prolapse of the anterior and / or posterior vaginal wall, additional stabilisation of the affected vaginal wall can be achieved with the **DynaMesh®-PRP visible** implants measuring 17 cm x 15 cm. The method here is the same with additional fixation of the enlarged middle mesh segment on the affected vaginal wall. The depth of the dissection in the vaginal wall towards the bladder neck depends on the type and nature of the prolapse and is determined by the surgeon intraoperatively. The implants can be cut accordingly to the correct dimensions.

To this end, the vagina should be taut in order to separate the bladder from the pelvic floor fascia on the anterior wall. The tension is very important in order to distinguish the tissue layers clearly. A Breisky retractor is recommended, for example, as this spreads out the tissue. The width should be tailored to the size of the defect.

The bladder should be dissected no further than the bladder neck. After fixing the PRP in the central area of the apex, the front part of the mesh can be placed on the exposed area. The length must be adapted accordingly. The caudal ends should be fixed in place with a monofilament (e.g., PDS). The all-important load-bearing fixation is positioned at the apex; the material placed here serves to reinforce the fascia and must therefore be crease-free.

In the case of dorsal defects (rectocele / enterocele), the septum is opened recto-vaginally and the rear vaginal wall is exposed. Here, too, the tissue should be pulled taut to achieve a vessel-free view under the vagina.

The extent of the dissection is based, in turn, on the defect and is made at the discretion of the surgeon. The **DynaMesh®-PRP visible** 17 cm x 15 cm is designed to enable the levator muscles to be reached. Again here the tissue must be adapted to the extent of the dissection. To avoid slipping, the distal points should be fixed in place with absorbable monofilament material.

Fixation

For fixation of the **DynaMesh®-PRP soft / visible** to the iliopectineal ligament and the cervical stump, non-absorbable suture material with the USP thickness 0 is recommended. Sharp needles (e.g., HS25) are recommended for use on the ligament and cervix. Absorbable, monofilament suture material with the USP thickness 2.0 should be used for fixation to the vaginal stump. We have no experience of using other fixation materials such as staplers or glues.

3. CONTRAINDICATIONS

All the usual contraindications relating to the use of implants for tissue support also apply to the application of **DynaMesh®-PRP soft / visible** mesh implants. **DynaMesh®-PRP soft / visible** mesh implants should not be used in patients with obvious infections, if the patient is currently or intends to become pregnant, or if the patient is still growing.

4. WARNINGS

DynaMesh®-PRP soft / visible mesh implants are intended for single use only. Re-use can lead to functional, immunological, toxicological or hygienic impairments. The user must be familiar with the surgical technique involved in the placement of mesh before using **DynaMesh®-PRP soft / visible**. The use of alloplastic mesh should only be undertaken after a careful consideration of classic surgical procedures has been made and the patient has been offered a comprehensive explanation.

If **DynaMesh®-PRP soft / visible** mesh implants are used in a contaminated wound, a subsequent infection is possible which may necessitate removal of the material. The known risks in prolapse surgery with the aid of alloplastic materials include erosions, infections, pain, bleeding, dyspareunia, organ perforation, voiding dysfunction, de novo urinary incontinence and shrinkage of the mesh.

The following complications may arise with meshes in general: seromas, haematomas, fistula formation, inflammations, recurrences or erosions. The patient should be advised that physical activities such as heavy lifting, jogging, riding or sexual intercourse should be resumed only after consultation with the doctor in attendance. Strictly sterile handling must be ensured up to implantation. The procedure must only be carried out in aseptic rooms. Damaged implants must no longer be used. Contact with sharp objects in particular must be avoided.

DynaMesh®-PRP soft / visible mesh implants must not be implanted if

- the expiry date has passed or
- the packaging has been damaged or
- there is reason to suspect the packaging is not intact.

The implant is supplied in a sterile double packaging. The inner packaging should not be opened until the mesh is about to be used. It is recommended that the mesh be handled only with sterile, unpowdered gloves and clean instruments. The surgeon is responsible for complications arising from an incorrect indication, an inadequate surgical technique or failure to ensure sterile conditions. In all cases, we recommend following up the procedure with appropriate post-operative treatment.

5. SIDE EFFECTS

There are no indications in the literature that the polyvinylidene fluoride used has any adverse side effects. As with all foreign bodies, **DynaMesh®-PRP soft / visible** mesh implants can have a negative effect on a pre-existing infection and can lead to temporary local irritation.

6. STERILITY

DynaMesh®-PRP soft / visible mesh implants are sterilised with ethylene oxide. **DynaMesh®-PRP soft / visible** must on no account be resterilised. The products must not be used if the packaging has been damaged or opened! Opened, unused products must be discarded!

7. STORAGE CONDITIONS

The recommended storage conditions are:

- Store in a dry place at room temperature,
 - Protect from moisture and the direct effect of heat,
 - Store only in the original packaging.
- Do not use after the expiry date!

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