Tailored Implants made of PVDF

DynaMesh[®]-ENDOLAP

- (de) Netzimplantat für die Hernienchirurgie aus PVDF-Monofilament
- (en) Mesh implant for hernia surgery, made of PVDF monofilament
- (ar) زرع شبكي خاص بجراحة الفتق من خيط وحيد مكون من ثنائي الفلوريد المتعدد الفيبيدي
- (cs) Síťkový implantát pro herniální chirurgii z PVDF-monofilamentu
- (da) Netimplantat af PVDF-monofilament til herniekirurgi
- (el) Εμφύτευμα πλέγματος για χειρουργική επέμβαση κήλης, κατασκευασμένο από μονονημάτιο PVDF
- (es) Implante de malla para cirugía de hernias fabricado en monofilamento de PVDF
- (fi) PVDF-monofilamentista valmistettu verkkoimplantti tyräkirurgiaan
- (fr) Filet prothétique pour la chirurgie des hernies, en monofil de PVDF
- (it) Reti per la chirurgia dell'ernia in monofilamento di PVDF
- (nl) Netimplantaat voor de herniachirurgie uit PVDF-monofilament
- (no) Nettingimplantat i PVDF-monofilament for brokkirurgi
- (pl) Siatka chirurgiczna z włókien monofilamentowych z polifluorku winylidenu (PVDF) do chirurgicznego leczenia przepuklin
- (pt) Implante de rede de monofilamentos de PVDF para a cirurgia de hérnias
- (ro) Implant de plasă din monofilament PVDF ptr.chirurgia herniei
- (ru) Сетчатый имплантат для герниопластики из ПВДФ мононитей
- (sk) Nettingimplantat i PVDF-monofilament for brokkirurgi
- (sv) Nätiimplantat för bräckkirurgi av PVDF-monofilament
- (tr) Herni cerrahisinde kullanılan PVDF Monofilament ağ implant
- (zh) 疝氣手術植入網片 (由 PVDF 單纖維製成)

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Instructions for Use

1. Description

DynaMesh®-ENDOLAP is a mesh implant for reinforcing connective tissue structures. **DynaMesh®-ENDOLAP** does not have any pharmacological effect of its own and is non-toxic. **DynaMesh®-ENDOLAP** is knitted using non-absorbable, biostable polyvinylidene fluoride monofilament. For handling and recognition purposes, the threads are partially coloured green or black. The coloured marking threads are used to improve placement of the mesh during the operation, in addition, the black threads enable postoperative recognition of the mesh in an MRI (Magnetic Resonance Imaging) scan.

The specially developed knitted structure and the optimal mesh size of the **DynaMesh®-ENDOLAP** are perfectly matched to the dynamometry of the groin region.

2. Fields of application

DynaMesh®-ENDOLAP serves to support the tissue and stabilise the fascial structures of the groin. **DynaMesh®-ENDOLAP** has been specially developed for the endoscopic (laparoscopic) repair of inguinal hernias using common minimally invasive surgical techniques (TEP and TAPP). **DynaMesh®-ENDOLAP** must not be used for the intraperitoneal onlay mesh technique (IPOM).

The mesh can be cut to individual requirements outside of the operation field using sharp scissors or a scalpel. The dimensions of the mesh should be such that it overlaps the hernia gap by approx. 3-5 cm on all sides.

It must be taken into consideration that, according to the current state of knowledge, 10 x 15 cm is the minimum recommended size required to cover the inguinal region. When using smaller meshes, there is an increased risk of recurrent hernias following endoscopic-laparoscopic treatment.

Fixation

The mesh can be fixed using either absorbable or non-absorbable suture material as well as absorbable or non-absorbable clips, staples or clamps. Adhesives can also be used. Based on literature and, as such, the current state of knowledge, the waiver of all types of fixation would seem possible, at least for total extraperitoneal hernioplasty (TEP). When using **DynaMesh®-ENDOLAP** as part of transabdominal preperitoneal hernioplasty (TAPP), fixation must be achieved in the above manner. The size of the clamps and the distances between them are to be selected by the surgeon based on empirical values.

3. Contraindications

All of the usual contraindications relating to use as an implant for tissue support also apply to the use of **DynaMesh®-ENDOLAP**. Among other things, these include:

- The mesh should not be used where there has been a previous enterotomy, due to the increased risk of infection.
- The mesh should be separated from the abdominal cavity by the peritoneum.
- The mesh is not suitable for use as a plug.
- Caution must be exercised with patients with obvious infections and pregnant women.
- A non-absorbable mesh should be not used for inguinal hernia repairs on growing patients and children.

4. Warnings

DynaMesh®-ENDOLAP is intended for single use only. Re-use can lead to functional, immunological, toxicological or hygienic impairments. The user must be familiar with the technique of inguinal hernia repair before using **DynaMesh®-ENDOLAP**. If **DynaMesh®-ENDOLAP** is used in a contaminated area, a subsequent infection is possible which may necessitate removal of the material. The implantation must be carried out under strictly sterile conditions. The application must be carried out only in aseptic rooms. Damaged implants may no longer be used. In particular, contact with sharp objects is to be avoided.

DynaMesh®-ENDOLAP must not be implanted if:

- the expiry date has been exceeded
- the packaging is damaged or
- there is reason to suspect the packaging is not intact.

The implant is supplied in sterile double packaging. The inner packaging should not be opened until immediately before the mesh is to be used. It is recommended that the mesh be handled only with sterile, unpowdered gloves (change gloves before use) and clean instruments.

The following complications may arise with meshes in general:

seromas, haematomas, recurrences and infections, as well as visceral adhesions, should the implant come into direct contact with the intestines. The surgeon is responsible for complications arising from an incorrect indication, an inadequate operating technique or failure to ensure sterile conditions. In all cases, it is recommended that the operation be followed up with appropriate post-operative treatment.

5. Side effects

There are no indications in literature that the polyvinylidene fluoride used has adverse side effects. As with all foreign bodies **DynaMesh®-ENDOLAP** may have a negative effect on a pre-existing infection and may lead to temporary local irritations.

6. Sterility

The product is sterilised with ethylene oxide. It is not permissible to resterilise **DynaMesh®-ENDOLAP**! The product must not be used if the packaging has been damaged or opened. Opened, unused products must be discarded!

7. Storage conditions

The recommended storage conditions are:

- Store in a dry place at room temperature.
- Protect from moisture and the direct effect of heat.
- Store only in the original packaging.
- Do not use after the expiry date.

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