
**DynaMesh<sup>®</sup> -IPST**  
**DynaMesh<sup>®</sup> -IPST visible**  
**DynaMesh<sup>®</sup> -IPST-D visible**  
**DynaMesh<sup>®</sup> -IPST-R**  
**DynaMesh<sup>®</sup> -IPST-R visible**

de	Netzimplantat	2
en	Mesh implant	3
ar	زرع شبكي	4
cs	Sítkový implantát	5
da	Netimplantat	6
el	Εμφυτεύσιμο πλέγμα	7
es	Implante de malla	8
fa	ایمپلنت	9
fi	Verkkoistuke	10
fr	Filet prothétique	11
he	שתל רשת	12
hu	Hálóimplantátum	13
it	Trapianto di rete	14
nl	Netimplantaat	15
no	Nettimplantat	16
pl	Implant siatkowy	17
pt	Implante têxtil	18
sk	Sieťkový implantát	19
sv	Nätimplantat	20
tr	Ağ implant	21
zh-cn	植入网片	22
zh-tw	植入網片	23

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**DynaMesh®-IPST /-IPST visible**  
**DynaMesh®-IPST-D visible**  
**DynaMesh®-IPST-R /-IPST-R visible**

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### DynaMesh®-IPST-R /-IPST-R visible

#### 1. DESCRIPTION

**DynaMesh®-IPST** is a mesh implant for reinforcing connective tissue structures with a penetration point for intestinal parts. **DynaMesh®-IPST** does not have any pharmacological effect of its own and is non-toxic. **DynaMesh®-IPST** is knitted using non-absorbable, biostable PVDF and polypropylene monofilament fibre. The parietal side is identified by filament ends which can also be used for fixation. Black-coloured thread material in the mesh structure enables post-operative recognition of the mesh in an MRI (Magnetic Resonance Imaging) scan.

#### 2. FIELDS OF APPLICATION

**DynaMesh®-IPST** is used in the prophylaxis and therapy of hernias and fascial defects in the abdominal wall after ostomy.

**DynaMesh®-IPST** is used to reinforce connective tissue structures and prevent prolapse of the diverted portion of the intestines.

It is used both in laparoscopy and in open surgery.

Common applications follow the intraperitoneal Onlay Mesh technique since the mesh is composed of different materials, which allows intraperitoneal application. Different surface structures (parietal side = polypropylene, visceral side = polyvinylidene fluoride) ensure ideal integration into the abdominal wall and reduced adhesion between the intestine and the mesh. Therefore, it is absolutely vital that the surfaces are aligned properly. The penetration point is provided with a funnel that must be directed toward the abdominal cavity. This ensures integration of the stoma loop in the funnel without increasing the risk of penetration. Furthermore, this position effectively prevents intestinal prolapse.

The diameter of the surgical mesh's flue must account for the size of the loop stoma. A tight fit is desirable to avoid a prolapsing of the stoma. An optimal application of the mesh is guaranteed when a visible widening of the flue is observed when guiding the loop stoma through it.

The parietal side is marked with green filament ends for clear identification. The mesh implant

must be placed such that these filament ends are oriented cranially and towards the abdominal wall when a **DynaMesh®-IPST-D** is used. For all other variants, the filament ends can be oriented cranially or caudally and towards the abdominal wall.

**DynaMesh®-IPST** is selected from the product range according to the application and the size of the defect.

If the variant **DynaMesh®-IPST-D** is used, attention should be paid to side specificity (left-sided/right-sided stoma), through which a particularly large overlap in a cranial as well as in a medial direction is ensured.

In the case of all **DynaMesh®-IPST**, the base area outside of the operating field and corresponding to individual requirements can be customised using sharp scissors or a scalpel.

Cutting towards the middle of the mesh should be avoided. If a parastomal hernia is being treated without releasing the diverted portion of bowel from the abdominal wall, the **DynaMesh®-IPST-R** variant with a prefabricated slit is available.

This slit should be closed in a non-absorbable manner after placement of the mesh. When cutting, the side identification by means of the green filament ends as well as adequate overlapping must remain ensured.

#### Fixation

For fixation of the mesh, non-absorbent suture material should be used. If fixation only involves the suture, a spacing of 1-2 cm must be maintained between the individual sutures, and at least 1 cm of mesh border should be added. Alternatively, a continuous suture technique can be used. The size of the suture depends on the surgeon's preference and the type of reconstruction but should not be less than the filament/suture size of 2-0.

Alternatively, clips, staplers or clamps can be used. However, care must be taken to ensure that fixation is performed not only in the peritoneum but also in the fascias or muscular structures. Experience has shown that spiral clips are most suitable for this.

#### 3. CONTRAINDICATIONS

All the usual contraindications relating to use as a tissue support implant also apply to the use of **DynaMesh®-IPST**.

These include:

- Due to the increased risk of infection, the mesh should not be used if the abdominal cavity has previously been contaminated with faeces.

- The mesh is not suitable for use as a plug. Caution must be exercised with growing children, obvious infections and pregnant women. **DynaMesh®-IPST** cannot be used on patients with an established allergy to polypropylene or PVDF (very rare).

#### 4. WARNINGS

**DynaMesh®-IPST** is intended for single use only. Re-use can lead to functional, immunological, toxicological or hygienic impairments.

The user must be familiar with the hernioplasty techniques applied before using **DynaMesh®-IPST**.

In the event of faulty alignment of the implant, strong adhesion to the intestine and infectious complications and even fistula formation are possible. Moreover, the mesh may fail to grow into the abdominal wall. If **DynaMesh®-IPST** is used in a contaminated wound, i.e. not in the abdominal cavity with all-round contact with the peritoneum, there is a possibility of subsequent infection, which may necessitate the removal of the material. Strictly sterile handling must be ensured up to implantation. The procedure must only be carried out in aseptic rooms. Damaged implants must no longer be used. Contact with sharp objects in particular must be avoided.

**DynaMesh®-IPST** must not be implanted if,

- the expiry date has passed or
- the packaging is damaged or
- there is reason to suspect the packaging is not intact.

The implant is supplied in a sterile double packaging. The inner sterile packaging should not be opened until the mesh is about to be used. It is recommended that the mesh be handled only with sterile, unpowdered gloves and clean instruments. The following complications may arise with meshes in general: seromas, haematomas, relapses, infections, visceral adhesions, mesh migration with fistula formation. The surgeon is responsible for complications arising from an incorrect indication, an inadequate surgical technique or failure to ensure sterile conditions. In all cases, we recommend following up the procedure with appropriate post-operative treatment.

#### 5. SIDE EFFECTS

There are no indications in the literature that the polypropylene or PVDF used has any adverse side effects. As with all foreign bodies, **DynaMesh®-IPST** may have a negative effect

on a pre-existing infection and may lead to temporary local irritations. Procedures involving intraperitoneal implantation include the general risk that points of adhesion will develop between the implant and the intestines. This applies, in particular, if the patient has undergone adhesiolysis prior to the implantation or if the peritoneum has been damaged in some other way.

#### 6. STERILITY

The product is sterilised with ethylene oxide.

**DynaMesh®-IPST** must on no account be resterilised!

The product must not be used if the packaging has been damaged or opened! Opened, unused products must be discarded!

#### 7. STORAGE CONDITIONS

The recommended storage conditions are:

- Store in a dry place at room temperature,
- Protect from moisture and the direct effect of heat,
- Protect against ultraviolet light,
- Store only in the original packaging,
- Do not use after the expiry date.

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